

CheckLite™ 250 Plus



Characteristics of kit

- CheckLite™ 250 Plus is a kit for microbial biomass assays based on the ATP measuring method.
- CheckLite™ contains thermostable firefly luciferase.
- CheckLite™ is characterized by its ease of use and rapid performance.
- CheckLite™ can detect 10,000 cells of coliforms/ml (1,000 cell/assay).

Composition of kit

(250 assays/kit)

Luciferin-luciferase reagent (lyophilized)	5
Reconstitution buffer (5.5 ml)	5
ATP releasing reagent (5.5 ml)	5

CheckLite™ 250 Plus should be stored at 2-8°C in a refrigerator. DO NOT FREEZE.

Manufacturer, Contact



Kikkoman Biochemifa Company

2-1-1, Nishi-shinbashi, Minato-ku, Tokyo 105-0003, Japan

Phone: +81-3-5521-5492 Fax: +81-3-5521-5498

E-mail: biochemifa@mail.kikkoman.co.jp

URL: <http://biochemifa.kikkoman.co.jp/>